

Head Office or Zone 1 Constantia Boulevard Constantia Kloof Ext 22, 1709

Postal Address: Private Bag X14, Weltevreden Park, 1715 Tel: 011 761 9000 Fax: 0866747666 Email: customerservice@sanbs.org.za www.sanbs.org.za

Antibody Identification and Titration Proficiency test Batch ABID07/19

Date: 2 August 2019

Antibody identification titration consensus score is obtained by using the mode of the submitted results. For batch ABID 07/19 there was two modes (1:4 and 1:128).

The titration values obtained during production was added to the results. Refer to FRM-QCL-115

Stability testing performed does not show any signs of deterioration of the antibody.

The mode calculated including the production titration results is 128, making the acceptable range for this batch 1:31 to 1:564

This was communicated and discussed with the PAC members.

Regards

Truscha Niekerk

Proficiency Testing Scheme Manager

Universal Blood Type:



ANTIBODY TITRE CONSENSUS





Set nr	1:1	1:2	1:4	1:8	1:16	1:32	1:64	1:128	1:256	1:512	1:1024	AINDUS
1				1							Regis	tration No. 2000/02/02/03
2				1								
3												
4												
5				1								
6					1							
7						1						
8								1				
9							1					
10												
11					1							
12												
13			1									
14				1								
15								1				
16								1				
17								1				
18								1				
19								1				
20							1					
21												
22							1					
23							1					
24							1					
25						1						
26												
27												
28												
29												
30												
31												
32												
33												
34												
35												
Total	0	0	1	4	2	2	5	6	0	0	0	0

Completed By:

Date: 25 07 1

Checked By: 4

Date: 02 08 2019